

## NOELLE<sup>®</sup> Birthing Simulator with birthing and resuscitation baby, with OMNI<sup>®</sup>.

### Product Information

Weight	35.99 kg
Brand	Gaumard
Item No.	W45111 [1012417]



### Description

Childbirth simulators have been an integral part of Maternal and Child health programs for more than 40 years. The childbirth simulator provides a comprehensive teaching system which combines the best of our patient care simulators with the advanced childbirth simulator. NOELLE<sup>®</sup> is designed to provide a complete birthing experience before, during and after delivery. The NOELLE<sup>®</sup> birthing simulator comes with four umbilical cords, four dilating cervices, two umbilical clamps, talcum powder, water soluble silicone lubricant, instruction manual and carrying bag.

The NOELLE<sup>®</sup> birthing simulator features:

- Full size articulating female
- Intubatable airway with chest rise
- IV arm for meds/fluids
- Removable stomach cover
- Practice Leopold Maneuvers
- Multiple fetal heart sounds
- Automatic birthing system
- Measure head descent and cervical dilation
- Multiple placenta locations
- Replaceable dilating cervices
- One articulating birthing baby with placenta
- One resuscitation baby with intubatable airway and umbilical catheter site

For a complete and realistic labor and birthing process use NOELLE<sup>®</sup> birthing simulator.

**Birthing**

**Birthing Birthing**

More product images

