

## CERTIFICATE OF COMPOSITION

### NON-STERIL LUBRICANT GEL

INGREDIENTS	CAS NUMBER	QUANTITY (w/v %)	FUNCTION
Deionized water	-	q.s 100	solvent
Hydroxy Ethyl Cellulose	9004-62-0	1.0-3.0	Gel Agent
Glycerin	56-81-5	10.0-25.0	Wetting agent
Mono Propylene Glycol	57-55-6	2,0-4,0	Wetting agent
Citric Acid	77-92-9	0.01-0.1	pH Adjuster
Methyl Hydroxybenzoate	99-76-3	0,05-0.2	Preservative

**Özlem ÇAĞRICI**  
Medikal AR-GE Yöneticisi  
Medical R&D Manager