

Product Specifications Sheet

Specification Sheet No	SSN-JAL007-005			
Product Name	Plain Swabs			
Product Description	Cotton Gauze Swabs			
Range	Premier			
Manufacturer	Allmed Medical Products Co., Ltd			
Material	100% Cotton Gauze, Warp and Weft threads			
Sterilisation	Sterile - ETO			
Issue Date	28 th April 2025	Version No	R02	

Intended use	Recommended for non-invasive procedures such as skin cleansing or wound debridement.
---------------------	--

Packaging						
Product code	Colour	Dimensions (cm)	Ply	Items per Pack	Qty- Inner	Qty – Shipper Carton
7000A	Blue	10 x 10	12	5	25	400 (16 PE bags of 25)
7001A	Blue	10 x 10	32	5	15	120 (8 PE bags of 15)
7001S	Blue	10 x 10	32	10	10	80 (8 PE bags of 10)
7015A	White	10 x 10	32	5	15	120 (8 PE bags of 15)
7015S	White	10 x 10	32	10	10	100 (10 PE bags of 10)
7099S	White	22.5 x 7.5	16	5	15	120 (8 PE bags of 15)
7170	White	5 x 5	8	5	30	180 (6 inner boxes of 30)
7172	White	5 x 5	12	5	30	180 (6 inner boxes of 30)
7174	White	7.5 x 7.5	12	5	30	180 (6 inner boxes of 30)
7176	White	10 x 10	8	5	30	180 (6 inner boxes of 30)
7178A	White	10 x 10	12	5	30	300 (10 PE bags of 30)
7178S	White	10 x 10	12	2	100	1000 (10 PE bags of 100)
7179	Green	7.5 x 7.5	12	5 (tied)	30	180 (6 inner boxes of 30)
7180	Green	10 x 10	12	5	25	150 (6 PE bags of 25)
7181	Blue	10 x 10	12	5	25	150 (6 inner boxes of 25)

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-02	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Device Directive 93/42/EEC		
Storage conditions	Store in a cool dry place. Keep away from heat and direct sunlight.		
Shelf life	5 years with the recommended storage conditions		
Contraindications/Warnings	None		
Markings	Premier Logo, Size, Quantity, Code Ref, Lot No, Expiry Date, CE Marking, Latex free Symbol, EO Sterile, Do not resterilise, Do Not Use If Package Is Damaged, Keep Dry Keep Away from Sunlight, Lot, Barcodes, Manufacturing Date, Expiry Date, EC-REP, Distributor name and address		
Additional information	The product is latex free.		
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	26/09/2024	Update to Mediq UK Brand replacing DSN533	Elizabeth Ogunkolade
R02	28/04/2025	Correct the Qty-Inner of 7001A & 7099S, qty-shipper carton of 7178A	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-02	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024